

A high concentration of the purest active and functional substances which, by leveraging the full potential of modern dermocosmetic technology, preserve the integrity of the encapsulated cosmetic until it is released.

**Dermakléb SkinCaps** come in daily doses and are stored in special dispenser containing **28 treatments** to take with you wherever you go

**Double Effect SkinCaps**  
Hyaluronic Acid in microspheres  
and Biopeptides  
**DAY TREATMENT**

**Eyes Lift SkinCaps**  
Vitamin A and Chrysin  
**NIGHT TREATMENT**

**Silver Mask**  
with Vitamin C  
and Silver Dag  
**NIGHT TREATMENT**  
WITH A PRACTICAL BRUSH



HOME  
PROGRAMME  
■■■■■  
**28**  
DAYS

- |   |   |   |  |
|---|---|---|--|
| 1 | REDUCED WRINKLES<br>right after the application     | 3 | ANTI-DARK CIRCLES<br>thanks to the oligopeptides           |
| 2 | +20%* OF COLLAGEN<br>and elastin, for a better tone | 4 | GLOWING SKIN<br>thanks to the Hyaluronic Acid microspheres |

\*Refers to the number of fibroblasts producing collagen and elastin All Dermakléb SkinCaps products are compatible with the Dermakléb product line.

# angel dental

implants • cosmetics • facial aesthetics

## Angel Dental

8-9 Queens Parade,  
Hove, BN3 8JG

T: 01273 019980

E: [info@theangeldentist.co.uk](mailto:info@theangeldentist.co.uk)

W: [www.theangeldentist.co.uk](http://www.theangeldentist.co.uk)

THE 3-IN-1 SYSTEM  
THAT REDEFINES THE RULES OF FACIAL TREATMENTS



# Aquaz

TECH LINE

AQUA3 is a **3-in-1** technology that combines the extraordinary detergent properties of **water-based peel containing sequential acids** with microdermabrasion and the **transdermal delivery** of cosmetic ingredients.

This is an extremely effective system that raises the bar for any facial treatment, taking the effective properties of active substances to a new level.



IT CLEANSSES YOUR FACE DEEPLY thanks to the water-based peel



IT REMOVES DEAD SKIN CELLS to obtain a more effective exfoliation



IT REBALANCES SEBUM AND MAKES COMEDONES DISAPPEAR with the power of vacuum



IT MOISTURISES YOUR SKIN DEEPLY thanks to the infuser handpiece



IT RESTRUCTURES YOUR TISSUES with transdermal delivery



**DIAMOND MICRODERMABRASION** it is effective on thickened skins, acne scars, deep wrinkles or stretch marks





Cleansing your skin only with cleansing milk and tonic is not enough to prepare it for treatment. The outermost stratum corneum, which thickens over the years, makes the skin less bright and soft and creates a barrier that prevents active substances to penetrate.

## AQUA3 is the new way of cleansing your skin.

### 1 WATER-BASED PEEL WITH SEQUENTIAL ACIDS

The aspirating handpiece, together with a water jet to which **Mandelic Acid** and **Salicylic Acid** have been added, removes superficial impurities even from thickened skins by aspirating dead skin cells, blackheads and whiteheads thanks to the powerful effect of vacuum.



### 2 DEEP MOISTURISATION WITH HYALURONIC ACID

By using the infuser handpiece with low-molecular weight **Hyaluronic Acid**, your skin is deeply moisturised and appears as glowing and smooth as never before.



### 3 RESTRUCTURING AND OXYGENATING YOUR TISSUES!

By making the most out of the transdermic delivery of the mesoporatic handpiece, your skin is restructured by the **Hyaluronic Acid** and Dermakléb serums.



**DEEP CLEANSING** thanks to the water-based peel [aspirating/infuser handpiece]



**EXTRA-EFFECTIVE TREATMENTS** with mesoporation



**MICRODERMABRASION** with Diamond grains [for scars and stretch marks too!]

The **diamond-grain handpiece** is a new device to treat the most critical areas, like the T-zone, in a selective and delicate manner, while at the same time reducing the appearance of scars and stretch marks by penetrating deep.

**9 different heads**, various levels of abrasive power, for an extra-exfoliating effect and amazing results.



### THE PROPERTIES OF DIAMOND TO RENOVATE EVERY TYPE OF SKIN

An ordinary water-based peel is not enough to treat largely thickened skins. With its unique diamond-grain handpiece, AQUA3 allows treating the most problematic skins by performing a controlled, non-invasive microdermabrasion that can remove the outermost stratum corneum, thus improving skin texture and restoring skin elasticity.

ALSO FOR STRETCH MARKS AND SCARS